

Digital Microscopes Comparison Chart

Digital Microscope Models	CAMERA			LENS		SOFTWARE		LIGHTING	
	Standard Camera	Camera Options	C-Mount Adapter	Lens Included	Lens Options	Software Included	Software Options	Light Included	Light Options
MiScope	640 x 480 (vga)			40x—140x		• Video ToolBox	• Video ToolBox PRO • Video Averager • Particle Analyzer	Built-in White LEDs	• UV • IR • PRO kit
MiScope-MP	1.3 MP			• 12-140x		• Video ToolBox PRO	• Video Averager • Particle Analyzer	Built-in White LEDs	• UV • IR • PRO kit
TetraCam	ZC-105 1.3 MP	ZC-305 INF X	ü	• 1-40x zoom lens		• Video ToolBox PRO	• FocusLock • Video Averager • Particle Analyzer	• Flourescent Lamp • LED flash-light	• LED ring light • Z-lite • Z-lite2
ZDM	ZC-105 1.3 MP	ZC-305 INF X	ü	• 9-65x zoom lens	• 15-105x zoom with detents • .5x re-ducer lens • 2x doub-ler lens	• Video ToolBox PRO	• FocusLock • Video Averager • Particle Analyzer	• LED Ring light • Portable: USB bar light	• Z-lite • Z-lite2 • Bottom Light Stand
ZDM 360	ZC-105 1.3 MP	ZC-305 INF X	ü	• 9-65 zoom lens with detents		• Video ToolBox PRO	• FocusLock • Video Averager • Particle Analyzer	• Z-lite2 OR • Dual Halo-gen Fiber Illuminator	
X-Loupe	Canon PowerShot SD 1000 7.1 MP			• 60x • 100x • 150x		• Canon camera software		Built-in White LEDs	