

1240 Sussex Turnpike, Suite 5, Randolph, NJ 07869 Phone: 973-933-2043 Fax: 973-933-2336 sales@zarbeco.com www.zarbeco.com

Microscope Accessories Comparison Chart

Product	Includes	Features	Included for Micro- scopes	Optional for Micro- scopes
Z-lite	Battery pow- ered (4AA batteries in- cluded)	 Two 20 lumen LED sources Goosenecks bendable to any position Portable Variable intensity light Compact 50,000 hour bulb life 		MiScopeMiScope-MPZDMZTC (TetraCam)
Z-lite2	AC powered	 Two 40 lumen LED sources Goosenecks bendable to any position Variable in- tensity light Compact 50,000 hour bulb life 	• ZDM 360 (choose Halo- gen Dual Fi- ber Illumina- tor)	MiScopeMiScope-MPZDMZTC (TetraCam)
LED Ring Light	AC powered	Variable intensity lightRing of LED lights	• ZDM	• ZTC (TetraCam)
Halogen Dual Fiber Optic Illumi- nator	AC powered	Variable intensity lightGoosenecks bendable to any position	• ZDM 360 (choose or Z- lite2)	
Flourescent Flood Lamp	AC powered	• 15 Watt Fluo- rescent fold- ing flood lamp	ZTC (TetraCam)	
ZDM Hard Case	• Foam insert	 Hard, water- tight, meets carry-on air- line require- ments 		• ZDM