

Software Comparison Chart

Product	Requires Video ToolBox, PRO, or Premier	Features	Microscopes included with	Cameras included with	Optional for Microscopes & Cameras	Operating System
Video ToolBox		<ul style="list-style-type: none"> • Capture still images and video • Labels, measurement,time stamp • Draw on live image 	<ul style="list-style-type: none"> • MISC • MISC-SW • MISC-IR/UV/UVIR 	<ul style="list-style-type: none"> • USB 2.0 VideoLink 		Windows XP/7/8
Video ToolBox PRO		<ul style="list-style-type: none"> • All Video ToolBox features • Advanced measurement tools • Calibration saved inside images • Export to EXCEL[®] • Transparent Overlays • Picture-in-Picture Comparison 	<ul style="list-style-type: none"> • MISC-MP2 • ZDM /ZDM-360/ ZTC • USB3 HDMI 	<ul style="list-style-type: none"> • ZC-203 • ZC-505 	<ul style="list-style-type: none"> • MISC • USB 2.0 VideoLink 	Windows XP/7/8
Video ToolBox Premier		<ul style="list-style-type: none"> • All Video ToolBox Pro features • Edge-detection assist • Precise cursor positioning • CAD overlays (optional plugin) • Export to CAD (optional plugin) • Side-by-side comparison • Opens large image files 			<ul style="list-style-type: none"> • All (except MISC) 	Windows XP/7/8
FocusLock	✓	<ul style="list-style-type: none"> • Live increased depth of focus enhancement 			<ul style="list-style-type: none"> • ZDM/ZDM-360/ZTC • ZC-203/ZC-505 	Windows XP/7/8
Video Aver-ager	✓	<ul style="list-style-type: none"> • Improves images in low light levels 			<ul style="list-style-type: none"> • All 	Windows XP/7/8
Particle Ana-lyzer	✓	<ul style="list-style-type: none"> • Counts particles/colors in im-ages 			<ul style="list-style-type: none"> • All 	Windows XP/7/8
Global Posi-tioning Sys-tem (GPS)	✓	<ul style="list-style-type: none"> • Includes GPS antenna to stamp position onto any still image or movie being captured 			<ul style="list-style-type: none"> • All 	Windows XP/7/8
Foot Switch	✓	<ul style="list-style-type: none"> • Foot switch and software for hands-free image capture 			<ul style="list-style-type: none"> • All 	Windows XP/7/8
MAC Software		<ul style="list-style-type: none"> • Capture still and moving images • Labels • Measurements • Time lapse video • Draw on live image 			<ul style="list-style-type: none"> • MISC/MISC-MP2 • ZDM/ZDM-360/ZTC • ZC-203/ZC-505 • USB3 HDMI 	MAC OS X